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**LIST OF CHEMICALS COMMONLY USED IN TISSUE
CULTURE LABORATORY**

Potassium nitrate(KNO ₃)	Glycine
Ammonium Nitrate (NH ₄ NO ₃)	L-Glutamine
Ammonium sulphate, (NH ₄) ₂ SO ₄	Casein hydrolysate
Magnesium sulphate (MgSO ₄ .7H ₂ O)	Myo-Inositol
Calcium chloride dihydrate (CaCl ₂ .2H ₂ O)	Indole-3-acetic acid (IAA), 97%
Potassium dihydrogen orthophosphate (KH ₂ PO ₄)	Naphthalene acetic acid (NAA), 99%
Potassium iodide (KI)	6-Benzylamino purine(BAP)
Boric acid (H ₃ BO ₃)	2,4-Dichlorophenoxy acetic acid (2,4-D)
Manganese sulphate (MnSO ₄ .4H ₂ O)	Kinetin 99%
Zinc sulphate (ZnSO ₄ .7H ₂ O)	Sucrose
Sodium molybdate (Na ₂ MoO ₄ .2H ₂ O)	Ethanol (absolute 100%)
Copper sulphate (CuSO ₄ .5H ₂ O)	Calcium hypochlorite
Cobalt chloride(CoCl ₂ .6H ₂ O)	Sodium hypochlorite
Aluminium chloride (AlCl ₃)	Sodium hydroxide
Ferrous sulphate (FeSO ₄ .7H ₂ O)	Hydrogen peroxide
Sodium EDTA (Na ₂ .EDTA.2H ₂ O)	Hydrochloric acid HCl
Polyvinyl pyrrolidone (PVP)	Silver nitrate
Nicotinic acid	Mercuric chloride
Pyridoxine hydrochloride HCl, 99.0-100%	Tween 20
Thiamine hydrochloride HCl, 98.5-100%	Agar
L-Ascorbic acid	Lysol

ANNEX - 2

**LIST OF COMMON EQUIPMENT AND GLASSWARE REQUIRED
FOR A PROFESSIONAL TISSUE CULTURE LABORATORY**

A. Media preparation equipment & glassware

Autoclave
Coarse balance for measuring in grams
Sensitive balance for measuring in milligrams
Distillation plant
Deep freezer
Refrigerator
Micro-wave oven
pH Meter
Hotplate with magnetic stirrer
Dispenser
Spatula
Metal/plastic racks
Basket
Wash bottle
Weighing bottle
Stirring bars and rod
Timer
Millipore filter system
Materials for closure of tubes or flasks
 Cotton plug
 Aluminum foil
 Metal caps plastic cap
Different types of glassware
 Measuring cylinder
 Volumetric flask
 Conical flask
 Bottle
 Screw neck reagent bottle
 Graduated bottle
 Petridishes with cover
 Test tube
 Pipette
 Pipette pump
 Pipette holder etc.

B. Culture isolation equipment

Laminar air flow cabinet
Sprit lamp! hot bead
Knives
Scissors
Forceps
Scalpel
Handle
Scalpel blades
Sterile petridish
Hard board

C. General equipment

Drying oven
Drying rack
Filter paper
Gloves
Tape (autoclavable)
Trolley
Detergent
Brushes

D. For culture room

Air condition
Shelves
Fluorescent tubes
Timer
Shaker

ANNEX - 3

COMPOSITION OF SELECTED PLANT TISSUE CULTURE MEDIA			
Compound	Concentration (mg/l)		
	White's	MS	B5
Inorganic			
NH ₄ NO ₃	--	1650	--
KNO ₃	80	1900	2527.5
MgSO ₄ .7H ₂ O	750	370	246.5
KH ₂ PO ₄	--	170	--
(NH ₄) ₂ SO ₄	--	--	134
Ca(NO ₃) ₂ .4H ₂ O	300	--	--
Na ₂ SO ₄	200	--	--
NaH ₂ PO ₄ .H ₂ O	19	--	150
KC	65	--	--
CaCl ₂ .2H ₂ O	--	440	150
KI	0.75	0.83	0.75
H ₃ BO ₃	1.5	6.2	3
MnSO ₄ .H ₂ O	5	22.3	--
MnSO ₄ .H ₂ O	--	--	10
ZnSO ₄ .7H ₂ O	3	8.6	2
Na ₂ MoO ₄ .2H ₂ O	--	0.25	0.25
MoO ₃	0.001	--	--
CuSO ₄ .5H ₂ O	0.01	0.025	0.025
CoCl ₂ .6H ₂ O	--	0.025	0.025
Fe ₂ (SO ₄) ₃	2.5	--	--
FeSO ₄ .7H ₂ O	--	27.8	--
Na ₂ .EDT A.2H ₂ O	--	37.3	--
Sequestrene 330 Fe ⁴	--	--	28

ANNEX - 4

Composition of commonly used Murashiage and Skoog (MS) medium,1962		
Components	Concentrations	
(Stock - 1) Macronutrients	mg/L	For 100 L
KNO ₃	1900	190.00 g
NH ₄ NO ₃	1650	165.00 g
KH ₂ PO ₄	170	17.00 g
MgSO ₄ , 7H ₂ O	440	44.00 g
(Stock-2)		
CaCl ₂ , 2H ₂ O	370	37.00 g
(Stock-3) Micronutrients		
FeSO ₄ , 7H ₂ O	27.8	2.78 g
Na ₂ -EDTA	37.3	3.73 g
(Stock-4)		
MnSO ₄ , 4H ₂ O	22.0	2.20 g
H ₃ BO ₃	6.2	0.62 g
ZnSO ₄ , 4H ₂ O	8.6	0.86 g
KI	0.83	0.083 g
Na ₂ MoO ₄ , 2H ₂ O	0.25	0.025 g
CuSO ₄ , 5H ₂ O	0.025	0.0025 g
CoCl ₂ , 6H ₂ O	0.025	0.0025 g
(Stock-5) Vitamins		
Glycine	2.00	0.200 g
Nicotinic acid	0.50	0.050 g
Pyridoxine-HCl	0.50	0.050 g
Thiamine-HCl	0.10	0.010 g
Myo-inositol	100.0	10.00 g

MERISTEM CULTURE TECHNIQUE**Meristem culture**

Meristem culture technique involves isolation and in vitro culture of shoot apices containing the apical dome with a few primordial leaves. In this method virus-free planting materials can be produced by micro section of the growing point, nourished in artificial media under controlled environment. This method is also known as apical-tip culture, shoot-tip culture or culture of shoot apices. Apical meristems are highly organized yet undifferentiated tissue of around 90-110 cells that are capable of rapid multiplication.

Advantages of meristem culture

Meristem culture is a specialized technique of plant tissue culture for generation of virus free plantlets for the following reasons:

- Generally healthy and vigorous shoot tips are selected resulting in rejuvenation of the plant.
- Virus particles within the plant system are translocated by the vascular system. Since the meristems are devoid of vascular connections virus particles can hardly reach meristem.
- Virus particles invading meristems by chance cannot sustain there because of the very high concentration of growth hormones and other biochemical substances in the meristematic dome.
- Meristemic zones are covered with folds of leaf primordia which protect them from infection by insect vectors visiting the plant.
- Plant growth results from very rapid division of meristems and generally any invading microorganism cannot make pace with this rate of multiplication leaving the meristem free from infections.

A PRACTICAL PROTOCOL FOR PRODUCTION OF POTATO PLANTLET

In vitro culture using apical meristem of potato for shoot proliferation and root induction in MS medium supplemented with different types and concentration of phytohormones can produce virus free potato plantlets. The technique used for primary establishment of meristem culture is described below.

Explant Collection and Preparation

- Grow potato plants preferably either in the growth chamber or other controlled conditions.
- Excise shoot tips of 25-35 days old potato plants with the help of a sharp blade; collect shoot tips in a bottle containing distilled water with few drops of Detol and few drops of Tween-20 (Polyoxyethylene-20), and quickly bring to the laboratory.
- Wash the explants 2-3 times with gradual change of sterile distilled water.
- Carry out surface sterilization, meristem isolation and inoculations, and other aseptic manipulations (including sub culture) in a running laminar airflow cabinet.

Switch on the cabinet half an hour before use and clean with 70% ethyl alcohol to reduce the chances of contamination. Cover all instruments (scalpels, needle, forceps, tiles, petri dishes etc.) with aluminum foil and sterilize using the steam sterilization method. While working, sterilize these instruments periodically using 70% ethyl alcohol dip and the flaming method inside the inoculation chamber. To ensure complete aseptic condition wipe both hands with 70% ethyl alcohol.

Surface Sterilization of Explants

- Within a running laminar airflow cabinet, transfer explants to a 250 ml sterilized conical flask.
- Surface sterilize the explants by dipping in 0.1% HgCl₂ solution with gentle shaking for 2-8 minutes followed by 3-5 times washing with sterile distilled water.

Isolation of Meristem

- After sterilization place the explants on the sterile tiles using sterile forceps. Hold the shoot tip in one hand under the stereo-microscope with the help of a pair of forceps and snap the immature leaves and leaf primordia with slight pressure from the needle.
- Then gently isolate the exposed meristem tip that appears as shiny dome with a sharp blade.

Inoculation Technique

- After deplugging of culture tubes, very carefully place a single excised meristem tip (0.3 mm) on the 'M' shape filter bridge of the culture tubes containing liquid MS medium supplemented with different kinds of plant growth regulators for the establishment of primary meristem.
- Flame the neck of the tubes with spirit lamp and then plug. After inoculation, label the culture tubes by glass marker.
- The culture tubes are ready for incubation.

ANNEX - 6


**FLOW CHART : MERISTEM CULTURE TECHNIQUE FOR
VIRUS-FREE PLANTLET PRODUCTION LEADING TO
TC-BASED SEED POTATO PRODUCTION**


Sl.	Step	Activity
1	Sprouting of tuber	Induce sprouting of breeder grade seed tubers
2	Growing of plant	Raise plants in growth cabinet/or in net-house
3	Collection of shoot tip	(a) Collect shoot tips and wash in distilled water containing few drops of Dettol and Tween 20 (b) Wash explants 2-3 times with distilled water
4	Disinfection of shoot tip <i>inside a Laminar Air Flow Cabinet</i>	(a) Surface sterilize with 0.1% HgCl ₂ for 5-8 minutes (b) Wash explants 4-5 times with sterile distilled water
5	Isolation of apical/ axillary meristem <i>inside a Laminar Flow Cabinet</i>	Expose the meristem, and gently isolate the exposed meristem (<i>which looks like a dome</i>) under a stereomicroscope
6	Inoculation <i>inside a Laminar Flow Cabinet</i>	Inoculate on agar MS medium solidified in test tube
7	Regeneration of plantlets	Incubate in culture room under fluorescent light
8	Virus indexing	Identify virus disease-free plantlets through ELISA test
9	Sub culturing	Multiplication of virus-free plantlets through repeated sub-culturing of single node cuttings
10	Hardening of plantlets	Hardening in TC-lab as well as net-house beds
11	Planting of plantlets in net-house beds at optimal spacing (25 × 10 cm)	The hardened plantlet or sprout-cuttings taken from the hardened plantlets can be planted in net-house beds for the production of pre-breeder seed potatoes (mini-tubers)

12	Multiplication of mini-tubers	Production of Breeder seed potatoes in net-house using pre-breeder seed potatoes (mini-tubers)
13	Production of Foundation seed potatoes	Production of Foundation seed potatoes in open field using normal seed potato production practices
	<i>Steps 4-6 are to be performed inside Laminar Flow Cabinet</i>	<i>Virus indexing through ELISA test is a must for identifying virus-free plantlets, leading to subsequent multiplication of only virus-free plantlets</i>
<i>Note : Step 3 should be started in Nov-Dec in order to get desired number of virus-free plantlets for hardening in net-house beds (step 10) during October-November next year.</i>		

ANNEX - 7

FLOW CHART : SPROUT CULTURE TECHNIQUE FOR VIRUS-FREE PLANTLET PRODUCTION LEADING TO TC-BASED SEED POTATO PRODUCTION

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